

AC INDUCTIVE PROXIMITY SENSORS



Two stainless steel medical syringes are shown, one slightly behind the other. Both syringes have a clear plastic plunger and a metal barrel with a CE marking. The syringes are connected to thin, clear plastic tubing.

P7C1B1805AO3B2	11	P7C2B1808AO3B2
P7C1B1805AC3B2	12	P7C2B1808AC3B2

5.8φ3x0.5	
Approx.161 g	Approx.160 g



CE

P7C1B1805AO3S2	11	P7C2B1808AO3S2
P7C1B1805AC3S2	12	P7C2B1808AC3S2

4φ 3x0.34	
Approx.156 g	Approx.146 g

P7C1B1805AOMB3	51	P7C2B1808AOMB3
P7C1B1805ACMB3	51	P7C2B1808ACMB3

Approx.69 g	Approx.68 g
-------------	-------------

